

Supporting Information

Evaluation and Optimization of High-Field Asymmetric Waveform Ion Mobility Spectrometry for Multiplexed Quantitative Site-specific N-glycoproteomics

Pan Fang,¹ Yanlong Ji,^{1, 3, 4} Ivan Silbern,^{1, 2} Rosa Viner,⁵ Thomas Oellerich,^{3, 4, 6} Kuan-Ting Pan,^{1, 3, 4, *} Henning Urlaub^{1, 2, *}

¹Bioanalytical Mass Spectrometry Group, Max Planck Institute for Biophysical Chemistry, 37077 Göttingen, Germany.

²Institute of Clinical Chemistry, University Medical Center Göttingen, 37075 Göttingen, Germany.

³Hematology/Oncology, Department of Medicine II, Johann Wolfgang Goethe University, 60590 Frankfurt am Main, Germany.

⁴Frankfurt Cancer Institute, Goethe University, 60596, Frankfurt am Main, Germany.

⁵Thermo Fisher Scientific, 95134 San Jose, CA, United States.

⁶German Cancer Consortium/German Cancer Research Center, 69120 Heidelberg, Germany

***For correspondence:**

kuan-ting.pan@kgu.de

henning.urlaub@mpibpc.mpg.de

Table of Content:

Supplementary methods

Supplementary Figures:

Figure S1. Effects of FAIMS CV settings on glycopeptide identification.

Figure S2. Distribution of the precursor charge state detected in an LC-MS/MS measurement using varying FAIMS CVs.

Figure S3. Influence of charge states and the number of sialic acid on the delta CV.

Figure S4. CV scans of synthetic peptides and glycopeptides from -100 to 0 V (in 1 V step).

Figure S5. Comparison of the best single CV for peptides without (blue) and with glycans (red).

Figure S6. Intensity distributions of glycopeptides bearing different glycan groups in various CVs.

Figure S7. Influence of glycan structure and composition on the FAIMS mobility of TMT-labeled DG75 N-glycopeptides.

Figure S8. Distributions of precursor mass, precursor m/z and glycan mass for identified TMT-labeled and unlabeled IgM glycopeptides across the tested CVs.

Figure S9. Distributions of precursor m/z, precursor mass, glycan mass, and peptide mass of glycopeptides identified in DG75 across tested CVs.

Figure S10. Effects of FAIMS on multiplexed quantitation of N-glycopeptides.

Figure S11. The distribution of signal-to-noise ratios of reporter ions across tested CV settings.

Supplementary Tables:

Table S1: List of all synthetic (glyco)peptides and their detected optimal CVs.

Table S2. Putative glycan types and branches.

Table S3. Summary of the optimized FAIMS CVs in this work and recent studies.

Supplementary methods:

Preparation of Unlabeled and TMT-labeled Peptides from IgM. Human IgM purified from human serum (Sigma) was reconstituted at the final concentration of 1 mg/ml in 0.05 M Tris-HCl, 0.2 M sodium chloride, 15 mM sodium azide, pH 8.0. RapiGest SF surfactant (Waters) was dissolved in 50 mM triethylammonium bicarbonate (TEAB) to a concentration of 1% (wt/wt). The prepared RapiGest solution was added to the protein solutions to a final concentration of 0.1% (wt/wt). The proteins were heated at 95 °C for 10 min. After cooling to room temperature (RT), the proteins were reduced and alkylated with 10 mM tris(2-carboxyethyl) phosphine hydrochloride (TCEP, 500 mM stock, Sigma) and 20 mM iodoacetamide at 37 °C for 60 min in the dark. The proteins were digested using trypsin at an enzyme to protein ratio of 1:50 at 37 °C overnight. The digested samples were acidified by adding TFA to a final concentration of 1% (vol/vol). The insoluble particles in the samples were removed by centrifugation at 14,000 × g for 10 min using a benchtop centrifuge. The resulting peptides were cleaned up with an Oasis HBL column and dried in a SpeedVac concentrator. The same digested peptides from IgM were aliquoted into two parts for unlabeled and TMT labeled samples. For TMT labeled IgM peptides, TMT labeling procedure (Thermo Scientific) was performed according to the manufacturer's instruction. After TMT labeling, the resulting peptides were cleaned up with an Oasis HBL column.

Cell Culture. DG75 cells (DSMZ no.: ACC 83) were cultured in RPMI medium (Invitrogen) supplemented with 10–20% (vol/vol) heat-inactivated fetal bovine serum (Invitrogen), penicillin/streptomycin (Invitrogen) and L-glutamine (Invitrogen) at 37 °C with 5% CO₂. Cells were grown to ~90% confluency before harvesting and provided as cell pellets aliquoted at about 1 × 10⁷ cells per tube.

Preparation of TMT-labeled Peptides from DG75 Cells. Cell pellets were lysed in the lysis buffer consisting of 4% SDS (w/v), 50 mM HEPES, pH 8.0. The samples were sonicated for 15 min (15 s on, 15 s off) using Bioruptor at 4 °C. After centrifugation at 14,000 × g for 15 min, the supernatants were collected, and the protein concentration were determined using Pierce BCA Protein Assay Kit (Thermo Scientific). The proteins were reduced and alkylated with 10 mM TCEP (500 mM stock, Sigma) and 20 mM iodoacetamide at 37 °C for 60 min in the dark. A mixture of magnetic beads, including Sera-Mag SpeedBeads with a hydrophilic surface (GE Healthcare, cat.no. 45152101010250, Magnetic Carboxylate Modified), and Sera-Mag SpeedBeads with a hydrophobic surface (GE Healthcare, cat.no. 65152105050250, Magnetic Carboxylate Modified) was rinsed twice with water on a magnetic rack prior to use. The beads were added to protein lysates at the optimal working ratio of 10:1 (wt/wt, beads to proteins). The required minimum bead concentration is 0.5 µg/µL in order to provide sufficient surface for the immobilization of aggregated proteins. We then added acetonitrile (ACN) to protein lysates to a final percentage of 70% (vol/vol). The samples were allowed to stay off the rack for 10 min at room temperature (RT), followed by resting on the magnetic rack for 2 min at RT. The supernatant was discarded, and the beads were then washed for three times with 80% (vol/vol) ethanol. Beads were resuspended in 50 mM TEAB containing sequencing grade modified trypsin (1:50 of enzyme to protein amount) and incubated at 37 °C for 4 h or overnight in a ThermoMixer with mixing at 800 rpm. After digestion, we placed the tubes on a magnetic rack for a few minutes and transferred the supernatant to a fresh tube. TMT labeling procedure (Thermo Scientific) was performed according to the manufacturer's instruction.

Enrichment of Glycopeptides from DG75 Samples using Zwitterionic Hydrophilic Interaction Liquid Chromatography (ZIC-HILIC). Digested peptides from DG75 were acidified by adding 10% (vol/vol) trifluoroacetic acid (TFA) to a final concentration of 1% followed by centrifugation at 14,000 × g for 20 min. The supernatant was dried in a SpeedVac concentrator. We re-dissolved the dried peptides in loading buffer consisting of 80% (vol/vol) ACN and 1% TFA and maintained the peptide concentration at around 4 mg/ml. Meanwhile, we weighted out the ZIC-HILIC beads (5 µm, Welch) according to the peptide to bead ratio of 1:50 (wt/wt) and washed three times using the loading buffer. We loaded all beads onto a 200 µL pipette tip pre-packed with coffee filter. Samples were loaded five times (at least 4 min each time) followed by three times wash with loading buffer. The retained glycopeptides were eluted with 100 µL 0.1% TFA twice. The collected eluates were dried in a SpeedVac concentrator.

IgM and Yeast Peptide Interference Model. Yeast protein prepared from *S. cerevisiae* cells were purchased from Promega (Cat. V7341). Yeast proteins in 6.5 M urea/50 mM Tris-HCl (pH 8) at a protein concentration of 10 mg/ml were thawed on ice. Protein reduction and alkylation were the same as mentioned above. After diluting urea to 1 M using 50 mM Tris-HCl (pH 8), trypsin (Promega) was added at a trypsin:protein ratio of 1:50 at 37 °C. After overnight incubation, the samples were acidified by adding TFA to a final concentration of 1% (vol/vol). The insoluble particles in the samples were removed by centrifuging at 14,000 × g for 10 min using a benchtop centrifuge. The samples were cleaned up with Oasis HBL columns. The resulted peptides were dried in a SpeedVac concentrator.

We separately labeled IgM digests with individual TMT6 reagents and pooled them together afterward with the ratio of 10:4:1:1:4:10. In contrast, yeast peptides were labeled with only the first three channels of TMT6 reagents (126, 127, and 128) and mixed equally (**Fig. 6A**). We then spiked the pooled IgM peptides into the yeast peptide mixture with an equal amount. The mixed samples were then cleaned up with Oasis HBL columns followed by analysis with either Glyco-SPS-MS3 or standard MS2 methods (see below).

LC-MS/MS Analysis. Peptides were resuspended in 5% ACN, 0.1% FA and subjected for LC-MS/MS analysis using Orbitrap Exploris 480 or Orbitrap Fusion Tribrid Mass Spectrometers (Thermo Fisher Scientific), both coupled to a Dionex UltiMate 3000 UHPLC system (Thermo Fisher Scientific). Unless noted otherwise, all samples were analyzed using a C18 trap column (3 cm long; inner diameter, 100 µm; outer diameter, 360 µm) and a home-made analytical column (ReproSil-Pur 120 C18-AQ, 1.9 µm pore size, 75 µm inner diameter, Dr. Maisch GmbH, 30 cm), the gradient at a flow rate of 300 nl/min, mobile phase A and B consisting of 0.1% (vol/vol) formic acid (FA) and 80% ACN, 0.08% FA, respectively. The one-hour gradient started at 10% B at 3 min, increased to 45% B at 47 min, and then to 90% B in 0.1 min. After washing with 90% B for 5 min, the column was re-equilibrated with 5% B. The two hours gradient started at 10% B at 3 min, increased to 45% B at 107 min, and then to 90% B in 0.1 min. After washing with 90% B for 5 min, the column was re-equilibrated with 5% B. The three hours gradient started at 10% B at 3 min, increased to 45% B at 166 min, and then to 90% B in 0.1 min. After washing with 90% B for 5 min, the column was re-equilibrated with 5% B.

The MS settings for the Glyco-SPS-MS3 method measured on Orbitrap Fusion are the same with previous publication¹. For standard MS2 method on Orbitrap Exploris 480, the following optimized MS settings were used. MS1 settings: Orbitrap Resolution-120 k, Scan Range (m/z)-350 - 2,000, Maximum injection time-50 ms, AGC target-300% (3e6), Microscans-1, RF Lens-40%, and Data Type-Profile. Filter MIPS-Peptide, Relax Restrictions when too few Precursor are Found-True. Filter Intensity Threshold: Minimum intensity-5000, Maximum Intensity-1E+20, Relative Intensity Threshold-20. Filter Charge State: Include Charge States-2-7, Include Charge State 1-False, Include Undetermined Charge States -False. Filter Dynamic Exclusion: Exclude after n Times-1, Exclude Isotopes-True, Exclusion Duration-22 s, Mass Tolerance Low-10 ppm, Mass Tolerance High- 10 ppm. Filter Precursor Selection Range: Mass Range-700-2000 m/z. Data Dependent Properties: Number of Dependent Scans- 20. MS2 settings: Isolation - window-1.6 m/z, Scan range mode-Auto normal, First mass-120, Normalized collision energy (%) -30/35/40, Detector type-Orbitrap, Orbitrap resolution-15 K, Maximum injection time-100 ms, AGC target-(50%)5e4, and Data type-Profile.

Data Analysis. For intact glycopeptide identification and quantification, .raw files were processed via GlycoBinder¹. Parameters used for pGlyco 2 include fully specific trypsin digestion with maximal two missed cleavage and mass tolerance for precursors and fragment ions of 10 and 20 ppm, respectively. We considered cysteine carbamidomethylation and TMT0 (or TMT6) on peptide N-termini and lysine residues as fixed modifications and methionine oxidation as a variable modification. The reviewed human protein database was downloaded from Swiss-Prot (March 2018, human, 20,303 entries). For the identification of IgM glycopeptides, we included only the sequences of human IGHM and IGJ in the FASTA file. Only GPSMs with PepScore ≥ 7 and GlyScore ≥ 8 reported by pGlyco 2 were used for following analysis. For the DG75 samples, we used the total FDR ≤ 2% for both the first and second database search.

For the IgM-Yeast interference model, we concatenated a reviewed *S.cerevisiae* protein database (Uniprot, 4,525 entries, March 2020, strain ATCC 204508 / S288c) with the sequences of human IGHM and IGJ protein for database search. FDR < 2% for both the first and second database search was used. A global analysis of all TMT ion ratios obtained for yeast and human peptides is shown for LC-MS/MS analyses performed with and without FAIMS. GPSMs of IgM glycopeptides that did not show reporter ion intensities in all TMT channels were filtered out.

Classification of Glycan Compositions. Glycan compositions were classified as (1) “oligomannose” if they consist of two HexNAc and more than five Hex; (2) “paucimannose” if they match the composition of Hex(1-4)Fuc(0-1)HexNAc(2), where the number of HexNAc is 2, the number of Fuc is 0 or 1, and the number of Hex is 1 to 4; (3) “short” if they match the composition of Hex(0-2)Fuc(0-1)HexNAc(1-3); (4) “complex” if they contain Hex-HexNAc antenna(e) attached to the common N-glycan core (Hex3HexNAc2); and (5) “hybrid” if they contain both Hex-HexNAc antenna(e) and extra Hex residues attached to the N-glycan core. There is no general composition for complex and hybrid types, as their structures are relatively more complex. We additionally include a “complex/hybrid” type for the glycan compositions that cannot be assigned as complex or hybrid types unambiguously.

We propose the potential number of branches of each glycan composition based on the glycan types. The number of branches in complex and hybrid type glycans was determined by the lowest possible number of antennae extended from the N-glycan core based on our understanding of human N-glycosylation biosynthesis. For the oligomannose type, the number of branches is 3. For the paucimannose and short glycan types, the number of branches is 1 or 2. Please see **Table S2** for the complete list of glycan compositions and their classifications.

Supplementary Figures:




Figure S1. Effects of FAIMS CV settings on glycopeptide identification. (A) Ratio of MS2 spectra containing 204.078 to total MS2 scans in a single LC-MS/MS run using varying CVs. Dashed lines represent the analyses without FAIMS. (B) Overlaps of unique TMT-labeled (blue) and unlabeled glycopeptides identified in LC-MS/MS runs using different CVs.




Figure S2. Distribution of the precursor charge state detected in an LC-MS/MS measurement using varying FAIMS CVs. (A) Charge state distribution of all MS1 features detected from TMT-labeled (left) or unlabeled (right) IgM samples across tested CVs. (B) Charge state distribution of identified glycopeptides with or without TMT-labeling across tested CVs.




Figure S3. Influence of charge states and the number of sialic acid on the delta CV. We detected differences (i.e., delta CVs) in the optimal CVs between TMT-labeled glycopeptides and their unlabeled counterparts. (A) The distribution of delta CVs for N-glycopeptides with different charge states. (B) The distribution of delta CV for N-glycopeptides with different numbers of sialic acids.




Figure S4. CV scans of synthetic peptides and glycopeptides from -100 to 0 V (in 1 V step). Ion signals of all (glyco)peptides with different charge states detected with varying CVs were extracted and shown in different colors. The m/z used for extracting signals are marked.




Figure S5. Comparison of the best single CV for peptides without (blue) and with glycans (red). We selected the (glyco)peptides sharing a common peptide sequence of ENISDPTSPLR for comparison. The best CV is the CV that provided the maximum intensity of each (glyco)peptide among all tested CVs. The glycan compositions are represented by the number of Hex, HexNAc, NeuAc, and Fuc. Each glycan composition was classified into different glycan types according to the criteria described in the Method section.




Figure S6. Intensity distributions of glycopeptides bearing different glycan groups in various CVs. To preclude the influence of peptide sequence, modifications and charge state, we considered only the triply charged IgM N-glycopeptides carrying the peptide sequence “YKNSDISSTR” (TMT6-labeled N-terminal and lysine residue on the TMT-labeled N-glycopeptides and no modification on the unlabeled glycopeptides). We then classified the identified N-glycopeptides into different groups based on their glycan type (see the classification of glycan compositions in supplementary methods). (A, E), number

of sialic acid moieties (B, F), number of possible branches (C, G), and number of glycan units (D, H). For each identified glycopeptide, the MS1 intensities detected with various CVs are normalized to their maximal intensity (100%) (i.e., normalize per row) and shown with a color scale. Four representative glycopeptides with different glycan types (dashed rectangles and arrows in Figure S2A and 2E) are also shown as intensity plots (I, J). Symbol nomenclature for glycans: Blue square: N-Acetyl-D-glucosamine, red triangle: fucose, green circle: mannose, yellow circle: galactose, purple diamond: sialic acid.



Figure S7. Influence of glycan structure and composition on the FAIMS mobility of TMT-labeled DG75 N-glycopeptides. We classified the selected DG75 N-glycopeptides into different groups based on their glycan type (see the classification of glycan compositions in supplementary methods), number of sialic acid moieties, number of possible branches, and number of glycan units. For each group (in columns), counts of unique N-glycopeptides identified across applied CVs (in rows) are shown. The color scale indicates the differences in counts as compared to the maximum number of identified unique N-glycopeptides in each group. The data include two measurement replicates.




Figure S8. Distributions of precursor mass, precursor m/z and glycan mass for identified TMT-labeled and unlabeled IgM glycopeptides across the tested CVs. All IgM glycopeptides are derived from five unique glyco-site-containing peptide sequences. Boxplots show the median (centerline), first and third quartiles (box edges) and $1.5 \times$ the interquartile range (whiskers).




Figure S9. Distributions of precursor m/z, precursor mass, glycan mass, and peptide mass of glycopeptides identified in DG75 across tested CVs. Boxplots show the median (centerline), first and third quartiles (box edges) and $1.5 \times$ the interquartile range (whiskers).




Figure S10. Effects of FAIMS on multiplexed quantitation of N-glycopeptides. The IgM-yeast mixture was analyzed by LC-FAIMS-MS2 or MS3 using varied CV settings. TMT reporter ratios were determined using the GlycoBinder pipeline. The ratio distributions of TMT channels with (red) and without (grey) yeast interference are showed in boxplots. Expected ratios are indicated. Boxplots show the median (centerline), first and third quartiles (box edges), and $1.5\times$ the interquartile range (whiskers). Outliers are not shown.




Figure S11. Distribution of signal-to-noise ratios of reporter ions across tested CV settings. Boxplots show the median (centerline), mean (squares), first and third quartiles (box edges) and $1.5 \times$ the interquartile range (whiskers). Outliers are not shown.

Supplementary Tables:

Table S1: List of all synthetic (glyco)peptides and their detected optimal CVs.

	Peptide sequences	Modification	Charge	Monoisotopic m/z	Peak CV (v)
P1	EVFVHPYSN N K	-	2+	610.31	-55
			3+	407.21	-44
P2	EVFVHPYS D K	-	2+	610.8009	-53
			3+	407.54	-42
GP1	EVFVHPYSN* K	GlcNAc	2+	711.8	-46
			3+	474.9	-47
GP2	EVFVHPYSN* K	5Hex4HexNAc2NeuAc	3+	1142.13	-45
			4+	856.85	-49
GP3	YGNVNETQNN* S FK	5Hex4HexNAc2NeuAc	3+	1240.49	-46
GP4	YGNVN* E TQNN S FK		2+	1860.73	-23
			3+	1240.49	-44
			4+	930.62	-70

Note: The amino acid difference between P1 and P2 is marked in red. The glycosylation sites are marked using a red star.

Table S2. Putative glycan types and branches.

Number of glycan units				glycan.Type	number.of.branches
Hex	HexNAc	Fucose	NeuAc		
0	1	1	0	short	1
0	2	0	0	short	1
0	2	1	0	short	1
1	2	0	0	paucimannose	1
1	2	1	0	paucimannose	1
2	2	0	0	paucimannose	1
2	2	1	0	paucimannose	1
3	2	0	0	paucimannose	2
3	2	1	0	paucimannose	2
4	2	0	0	paucimannose	2

4	2	1	0	paucimannose	2
5	2	0	0	oligo-mannose	3
5	2	1	0	oligo-mannose	3
6	2	0	0	oligo-mannose	3
7	2	0	0	oligo-mannose	3
8	2	0	0	oligo-mannose	3
8	2	1	0	oligo-mannose	3
9	2	0	0	oligo-mannose	3
9	2	1	0	oligo-mannose	3
10	2	0	0	oligo-mannose	3
11	2	0	0	oligo-mannose	3
12	2	0	0	oligo-mannose	3
1	3	0	0	short	1
1	3	1	0	short	1
2	3	0	0	short	2
2	3	1	0	short	2
3	3	0	0	hybrid	2
3	3	0	1	hybrid	2
3	3	1	0	hybrid	2
3	3	1	1	hybrid	2
3	3	2	0	hybrid	2
4	3	0	0	hybrid	2
4	3	0	1	hybrid	2
4	3	1	0	hybrid	2
4	3	1	1	hybrid	2
4	3	1	2	hybrid	2
4	3	2	0	hybrid	2
4	3	2	2	hybrid	2
4	3	3	0	hybrid	2
5	3	0	0	hybrid	3
5	3	0	1	hybrid	3
5	3	1	0	hybrid	3
5	3	1	1	hybrid	3
5	3	2	0	hybrid	3
5	3	2	2	hybrid	3
5	3	3	0	hybrid	3
6	3	0	0	hybrid	3
6	3	0	1	hybrid	3
6	3	0	2	hybrid	3
6	3	1	0	hybrid	3
6	3	1	1	hybrid	3
6	3	1	2	hybrid	3
6	3	2	1	hybrid	3
6	3	3	0	hybrid	3
7	3	0	0	hybrid	3
8	3	0	0	hybrid	3
2	4	0	0	hybrid/complex	2
2	4	1	0	hybrid/complex	2
3	4	0	0	hybrid/complex	2/3
3	4	0	1	hybrid/complex	2/3
3	4	1	0	hybrid/complex	2/3
3	4	1	1	complex	2
3	4	3	0	complex	2
4	4	0	0	hybrid/complex	2/3
4	4	0	1	hybrid/complex	2/3
4	4	0	2	hybrid/complex	2/3
4	4	1	0	hybrid/complex	2/3
4	4	1	1	hybrid/complex	2/3
4	4	2	0	hybrid/complex	2/3
4	4	2	1	hybrid/complex	2/3
4	4	3	0	hybrid/complex	2/3

4	4	3	1	hybrid/complex	2/3
4	4	4	0	hybrid/complex	2/3
5	4	0	0	hybrid/complex	2/4
5	4	0	1	complex	2
5	4	0	2	complex	2
5	4	0	3	complex	2
5	4	1	0	hybrid/complex	2/4
5	4	1	1	complex	2
5	4	1	2	complex	2
5	4	2	0	hybrid/complex	2/4
5	4	2	1	complex	2
5	4	2	2	complex	2
5	4	3	0	hybrid/complex	2/4
5	4	3	1	complex	2
5	4	4	0	complex	2
5	4	5	0	complex	2
6	4	0	0	hybrid/complex	2/4
6	4	0	1	hybrid/complex	2/4
6	4	1	0	hybrid/complex	2/4
6	4	1	1	hybrid/complex	2/4
6	4	2	0	hybrid/complex	2/4
6	4	3	0	hybrid/complex	2/4
7	4	0	0	hybrid/complex	2/4
7	4	0	1	hybrid/complex	2/4
7	4	1	0	hybrid/complex	2/4
7	4	1	1	hybrid/complex	2/4
7	4	2	0	hybrid/complex	2/4
3	5	0	0	complex	2/3
3	5	1	0	complex	2/3
3	5	3	1	complex	2/3
3	5	4	0	complex	2/3
4	5	0	0	complex	2/3
4	5	0	1	complex	2/3
4	5	0	2	complex	2/3
4	5	1	0	complex	2/3
4	5	1	1	complex	2/3
4	5	2	0	complex	2/3
4	5	3	0	complex	2/3
4	5	4	0	complex	2/3
5	5	0	0	complex	2/3
5	5	0	1	complex	2/3
5	5	0	2	complex	2/3
5	5	1	0	complex	2/3
5	5	1	1	complex	2/3
5	5	1	2	complex	2/3
5	5	2	0	complex	2/3
5	5	2	1	complex	2/3
5	5	3	0	complex	2/3
5	5	3	1	complex	2/3
5	5	5	0	complex	2/3
6	5	0	0	complex	2/3
6	5	0	1	complex	2/3
6	5	0	2	complex	2/3
6	5	1	0	complex	2/3
6	5	1	1	complex	2/3
6	5	1	2	complex	2/3
6	5	2	0	complex	2/3
6	5	2	1	complex	2/3
6	5	3	0	complex	2/3
6	5	4	3	complex	2/3
7	5	1	0	complex	2/3

3	6	2	0	complex	2/4
3	6	3	0	complex	2/4
3	6	4	0	complex	2/4
4	6	0	0	complex	4
4	6	0	1	complex	4
4	6	1	0	complex	4
5	6	0	0	complex	3/4
5	6	0	1	complex	3/4
5	6	0	2	complex	3/4
5	6	1	0	complex	3/4
5	6	1	1	complex	3/4
5	6	2	0	complex	3/4
5	6	2	1	complex	3/4
5	6	3	1	complex	3/4
5	6	5	0	complex	3/4
6	6	0	0	complex	3/4
6	6	0	1	complex	3/4
6	6	2	2	complex	3/4
6	6	3	1	complex	3/4
6	6	4	1	complex	3/4
7	6	4	0	complex	2/3/4
9	6	0	1	complex	3/4
10	6	1	1	complex	4
10	6	2	1	complex	4
3	7	2	0	complex	4/5
3	7	4	0	complex	4/5
4	7	1	1	complex	4/5
4	7	3	0	complex	4/5
5	7	1	1	complex	4/5
7	7	2	1	complex	5
8	7	1	2	complex	5
10	7	0	0	complex	5
11	7	1	0	complex	5
3	8	0	0	complex	5
10	8	3	1	complex	5

Table S3. Summary of the optimized FAIMS CVs in this work and recent studies

	Objects	Samples	Labeling	Peak CV	Double CVs	Triple CVs	Instrument	Ref.
1	Unmodified peptides	K562	unlabel	-60	-50/-70	-50/-65/-85	Orbitrap Fusion Lumos Tribrid; FAIMS pro	Anal. Chem. 2018, 90, 9529–9537.
		Rat brain/yeast/E.coli			-	-		
2	Unmodified peptides	Yeast cells	TMT6 label	-50	-50/-70	-40/-60/-80		Anal. Chem. 2019, 91, 4010–4016
		Saccharomyces cerevisiae, BY4742(TKO standard)	unlabel	-50	-	-		
3	O-linked glycopeptides	Burkholderia species (Burkholderia cenocepacia, Burkholderia Dolosa, and Burkholderia ubonensis),	Unlabel, Hex-HexNAc-HexNAcSuc-Hex-HexNAc-HexNAc	-30	-	-25/-35/-45		DOI: 10.1021/acs.jproteome.0c00565
	N-linked glycopeptides	Campylobacter fetus subsp. Fetus	Unlabel, Hex-HexNAc4-diNAcBacHexNAc5-diNAcBac	-20	-	-		
4	Phosphopeptides	phospho-enriched HCT116 lysate (human)	TMTpro/TMT10	-80	-70/-80	-40/-60/-80		J. Proteome Res. 2020, 19, 554–560
	Phosphopeptides	phospho-enriched tryptic rat brain	unlabel	-50	-45/-65	-40 -55 -70 V		Anal. Chem. 2020, 92, 24, 15959–15967
5	Cross-linked peptides	Eight Protein Mixture, HEK293T cell lysate	DSSO, unlabel	-60	-50/-65; -50/-70	-40/-55/-70; -50/-60/-75		Anal. Chem. 2020, 92, 10495–10503
			DSS, unlabel	-50	-40/-55; -50/-60	-40/-55/-65		
6	N-Linked glycopeptides	IgM; glycopeptides enriched in DG75	TMT6	-45; -50	-45/-55; -40/-50	-40/-50/-60	Exploris 480; FAIMS pro	This study
			unlabel	-45; -65	-45/-65	-45/-50/-65		

Supplementary references:

1. Fang, P.; Ji, Y.; Silbern, I.; Doebele, C.; Ninov, M.; Lenz, C.; Oellerich, T.; Pan, K. T.; Urlaub, H. *Nat. Commun.* **2020**, *11*, 5268.