

Table 1S. Local standard troponin tests and cutoffs for MI diagnosis

<b>Center</b>	<b>Troponin test</b>	<b>MI Cut-Off</b>
Heidelberg	Hs TnT, Elecsys, Roche Diagnostic	14 ng/l
UKB, Berlin	Hs TnT, Elecsys, Roche Diagnostic	14 ng/l
CVK, Berlin	AQT-Test POCT, Radiometer hsTnT, Elecsys, Roche Diagnostic	30 ng/l 50 ng/l
Frankfurt	Hs TnT, Elecsys, Roche Diagnostic	14 ng/l
Bad-Nauheim	Hs TnT, Elecsys, Roche Diagnostic	14 ng /l
Mayen	TnI Ortho Clinical Diagnostics and from 19.4.16 TnI, LOCI, Siemens	50 ng/l
Wien	TnI, LOCI, Siemens	45 ng/l
Calais	TnI, Access, Beckman and Coulter	30 ng/l (97.5th %le)
Vilnius	Hs TnI, Architect, Abbott	for men 34,2 ng/l for women 15,6 ng/l
Budapest	Hs TnT, Elecsys, Roche Diagnostic (Cobas e411)	14 ng/l
High Wycombe	Hs TnI, Architect, Abbott	for men 34,2 ng/l for women 15,6 ng/l
Zollichberg, Zurich	TnI-Ultra, Centaur, Siemens	40 ng/l
Aarau	TnI, LOCI, Siemens	45 ng/l
Berlin Hedwigshöhe	Hs TnI, Architect, Abbott	15 ng/l
Dijon	TnI, Vista, Siemens	100 ng/l
Ankara	TnI, Access, Beckman and Coulter	40 ng/l (99th %le)
Tulln	Hs TnT, Elecsys, Roche Diagnostic	14 ng/l
Montpellier	Hs TnT, Elecsys, Roche Diagnostic (Cobas 8000/e602 analyzer)	14 ng/l